

ABOUT SSP AND AUS SYSTEMS

Our commitment to serving not just as a product developer and manufacturer, but as the partner our customers require to help them solve their specialty design challenges is what shapes the focus of Specialty Surgical Products, Inc. (SSP). It has always been our goal to offer the highest quality surgical devices, such as silicone silo bags and pediatric anal dilators with improved design and effectiveness to the pediatric surgery community. We hold to these same exceptional standards for all our products.

With this in mind, SSP has partnered with Aus systems to offer you the next generation suction rectal biopsy system for the diagnosis of Hirschsprung's Disease. For over 12 years Aus systems has manufactured and distributed the Aus biopsy system—a reusable suction rectal biopsy instrument for the diagnosis of Hirschsprung's Disease.

Although this instrument has been used successfully during this time, it has required ongoing maintenance and servicing by the manufacturer to keep the blades sharp and the suction consistent to ensure successful biopsies. This often requires a 4 – 6 week turn around resulting in downtime and ongoing costs to the hospital.

Aus systems has now developed the next generation suction rectal biopsy system that combines a single-use capsule procedure pack and a reusable handpiece. Capsules are supplied sterile and fully assembled with an ultra-sharp blade and airtight seal. Operators can now expect precision cutting and suction consistency for each biopsy.

Specialty Surgical Products, Inc. as the exclusive U.S. distributor of rbi2, is proud to offer you another fine product for your practice. Call any of us toll free, 888.878.0811 to discuss our products in more detail, discuss a product design you'd like to have developed, or to place your orders.

Sincerely yours,

Sherry J. Null
VP Marketing & Sales



SPECIALTY SURGICAL PRODUCTS, INC.

Toll Free 888.878.0811

406.961.0102 • Fax: 406.961.0103 • 1131 North U.S. Hwy 93 • Victor, MT 59875

Email: info@ssp-inc.com, Web: www.ssp-inc.com

M209A 09/06