


SILICONE PERCUTANEOUS BILIARY DRAINAGE CATHETER



SPECIFIC DESIGN POINTS:


- **Soft silicone elastomer** properties enhance patient comfort.
- **Radiopaque marker** at distal end and catheter tip for improved x-ray visualization and positioning.
- **Catheter sizes** from 8 to 24 French.
- **Multiple holes** for efficient drainage.
- **Luer lock adapter** to facilitate drainage connection.
- **Latex free.**

SILICONE PERCUTANEOUS BILIARY DRAINAGE CATHETER

| CATALOG # | DESCRIPTION | SIZE | DIMENSIONS* |
|-----------|---|-------|---------------------|
| SBC08-1 | Silicone Biliary Drainage Catheter, 8 FR | 8 FR | ID .070" : OD .105" |
| SBC10-1 | Silicone Biliary Drainage Catheter, 10 FR | 10 FR | ID .087" : OD .131" |
| SBC12-1 | Silicone Biliary Drainage Catheter, 12 FR | 12 FR | ID .105" : OD .158" |
| SBC14-1 | Silicone Biliary Drainage Catheter, 14 FR | 14 FR | ID .122" : OD .184" |
| SBC16-1 | Silicone Biliary Drainage Catheter, 16 FR | 16 FR | ID .140" : OD .210" |
| SBC18-1 | Silicone Biliary Drainage Catheter, 18 FR | 18 FR | ID .157" : OD .236" |
| SBC20-1 | Silicone Biliary Drainage Catheter, 20 FR | 20 FR | ID .175" : OD .263" |
| SBC22-1 | Silicone Biliary Drainage Catheter, 22 FR | 22 FR | ID .192" : OD .288" |
| SBC24-1 | Silicone Biliary Drainage Catheter, 24 FR | 24 FR | ID .210" : OD .315" |


All catheters are supplied 24" (61cm) lengths.

*Dimensions have the following tolerances:
 8-22FR +/- .005" ID; +/- .006" OD
 24FR +/- .003" ID; +/- .002" OD



Order for Dr.: _____
 Patient: _____
 Surgery Date: _____

| Catalog #: | Qty: |
|------------|-------|
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |



◆ Luer lock connector supplied sterile with each catheter.

Product supplied sterile. Each device supplied with Luer Lock Connector.

Limited Warranty, Limitation of Liability and Disclaimer of Other Warranties: SSP warrants that reasonable care was used in the manufacture and production of this product. Because SSP has no control over the conditions of use, patient selection, surgical procedure, post-surgical stresses, or handling of the device after it leaves our possession, SSP does not warrant either a good effect or against an ill effect following its use. SSP shall not be responsible for any incidental or consequential loss, damage or expenses directly or indirectly arising from the use of this product. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to any implied warranties of merchantability or fitness for use.



Toll Free 888.878.0811

406.961.0102 · Fax: 406.961.0103 · 1131 North U.S. Hwy 93 · Victor, MT 59875
 email: info@ssp-inc.com