

Dermaspan™ Low Pole Plus Expander



Winged Needle Infusion Set conveniently supplied sterile with each expander for incremental filling with your delivery system of choice.

Designed by a surgeon for surgeons


Designed in collaboration with Dr. Dennis Hammond, the Dermaspan™ breast tissue expander is ideal for breast reconstructive procedures.

- ▼ Suturable, durable tabs prevent migration and rotation
- ▼ More consistent and aggressive texture for optimal shell-tissue interface
- ▼ Smooth surface with tabs available
- ▼ Low pole design places apex or expansion where needed the most
- ▼ Designed for non-interference with inframammary fold development
- ▼ Magnetic Precision Point™ easily locates the injection port
- ▼ No palpable rings or ridges for improved patient comfort during expansion
- ▼ Reinforced base facilitates placement and expander “roll out”



SPECIALTY SURGICAL PRODUCTS, INC.

Dermaspan™ Low Pole Expanders



Order for Doctor: _____

Patient: _____

Surgery Date: _____

| Catalog No: | Quantity: |
|-------------|-----------|
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |
| _____ | _____ |

Product supplied sterile. Each device supplied with a Precision Point™ locator and Winged Needle Infusion Set. Please see package insert for additional product information.

Custom shapes and sizes also available

Dermaspan and Precision Point are trademarks of Specialty Surgical Products, Inc.

| LOW POLE PLUS, TEXTURED SURFACE WITH LATERAL EXTENSIONS | | | | |
|---------------------------------------------------------|-----------------------|---------|----------|--------------|
| Catalog No. | Suggested Fill Volume | A Width | B Length | C Projection |
| Mid-Height | | | | |
| LPP11 | 250 - 300 | 11.0 | 9.0 | 5.0 - 5.7 |
| LPP12 | 310 - 370 | 12.0 | 9.8 | 5.2 - 6.0 |
| LPP13 | 380 - 450 | 13.0 | 10.6 | 5.5 - 6.3 |
| LPP14 | 460 - 550 | 14.0 | 11.2 | 6.0 - 7.3 |
| LPP15 | 500 - 600 | 15.0 | 11.5 | 6.4 - 7.9 |
| LPP16 | 580 - 690 | 16.0 | 12.0 | 6.6 - 8.2 |
| Full-Height | | | | |
| LPP-FH11 | 250 - 300 | 11.0 | 10.2 | 5.0 - 5.7 |
| LPP-FH12 | 310 - 370 | 12.0 | 11.0 | 5.2 - 6.0 |
| LPP-FH13 | 360 - 430 | 13.0 | 11.9 | 5.4 - 6.3 |
| LPP-FH14 | 480 - 575 | 14.0 | 12.9 | 6.1 - 7.3 |
| LPP-FH15 | 600 - 720 | 15.0 | 13.8 | 6.4 - 7.9 |
| LPP-FH16 | 750 - 900 | 16.0 | 14.7 | 6.9 - 8.5 |

Volume in cc, Dimensions in cm

| ACCESSORIES | |
|-------------|---------------------------------------------------------------------|
| Catalog No. | Description |
| DET01 | Dermaspan Expander Templates, non-sterile, set (request a FREE set) |
| PPL01 | Precision Point™ Locator, sterile, each |
| WNIS01 | Winged Needle Infusion Set, sterile, set of 5 each |

- ▼ **NOMINAL DIMENSIONS**
Except as noted, expander dimensions are measured at the minimum suggested fill volume with expander placed on a flat surface. Variations may occur, and not every expander will conform exactly to the dimensions given.
- ▼ **MAGNETIC STRENGTH**
This device contains a strong rare-earth magnet that has been tested to be detectable through 45mm phantom tissue. This magnet is contraindicated in MRI procedures and where the magnetic field may affect other polarized devices (e.g., drug infusion devices, pacemakers, artificial sensing devices and similar type products).
- ▼ **LIMITED WARRANTY, LIMITATION OF LIABILITY AND DISCLAIMER OF OTHER WARRANTIES**
SSP warrants that reasonable care was used in the manufacture and production of this product. Because SSP has no control over the conditions of use, patient selection, surgical procedure, post-surgical stresses, or handling of the device after it leaves our possession, SSP does not warrant either a good effect or against an ill effect following its use. SSP shall not be responsible for any incidental or consequential loss, damage or expenses directly or indirectly arising from the use of this product. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including but not limited to any implied warranties of merchantability or fitness for use.



SPECIALTY SURGICAL PRODUCTS, INC.

1131 NORTH US HIGHWAY 93 • VICTOR, MT 59875 • 888.878.0811 • 406.961.0102 • 406.961.0103 FAX

WWW.SSP-INC.COM

M0651 09/11