

Lipofoam™ Fine Cell Foam Sheeting



Non-gel

| Catalog No. | Description |
|-------------|------------------------|
| LPF01 | Lipofoam Foam Sheeting |

Lipofoam is a trademark of Specialty Surgical Products, Inc.

Lipofoam™ is fabricated from the highest quality fine cell polyurethane foam for superior performance and patient comfort. It is designed to give uniform smooth compression, and may be inserted between the post-operative compression garment and the suctioned area, or adhered to skin surface with the use of appropriate adhesive spray. Lipofoam is easily cut to any desired size or shape, making it ideal for arms, chins, abdomen, thighs or any area requiring uniform healing.

- ▼ Finest quality medical grade polyurethane foam for patient comfort and consistent compression distribution
- ▼ Gel-free for greater flexibility in placement and patient movement
- ▼ Latex-free
- ▼ Convenient packaging of ten (10) size 8 x 11-inch sheets per box, packaged in two (2) clean pouches of five (5) each



SPECIALTY SURGICAL PRODUCTS, INC.

1131 NORTH US HIGHWAY 93 • VICTOR, MT 59875 • 888.878.0811 • 406.961.0102 • 406.961.0103 FAX

WWW.SSP-INC.COM

M120E 09/11