








SPECIALTY SURGICAL PRODUCTS, INC.

Softspan™ Extremity Expander Wall Chart

Page 1 of 2

Expander Shape	Remote Port Cat #		Integrated Port Cat #		Suggested Fill Volume	A Length	B Width	C Projection
	Smooth Shell	Textured Shell	Smooth Shell	Textured Shell				
	SRE-FO-025S*	SRE-FO-025*	n/a	SIE-FO-65	25-62cc	6.5cm	5.0cm	1.5-4.0cm
	SRE-FO-035S*	SRE-FO-035*		SIE-FO-76	35-89cc	7.5cm	6.0cm	1.5-4.3cm
	SRE-FO-060S	SRE-FO-060		SIE-FO-86	60-150cc	8.5cm	6.5cm	1.5-5.2cm
	SRE-FO-075S	SRE-FO-075		n/a	75-188cc	10.0cm	8.0cm	1.5-5.5cm
	SRE-FK-025S	SRE-FK-025	n/a	SIE-FK-74	25-62cc	7.0cm	4.5cm	1.5-3.5cm
	SRE-FK-080S	SRE-FK-080		SIE-FK-106	80-200cc	10.0cm	6.0cm	2.0-5.5cm
	SRE-FK-140S	SRE-FK-140		SIE-FK-127	140-350cc	12.5cm	7.5cm	2.5-6.5cm
	SRE-FK-270S	SRE-FK-270		n/a	270-675cc	14.0cm	9.0cm	4.0-8.0cm
	n/a	n/a	SIE-LR-064S	SIE-LR-064	60-120cc	6.0cm	4.0cm	3.5-5.0cm
	SRE-LR-160S	SRE-LR-160	SIE-LR-106S	SIE-LR-106	160-320cc	10.0cm	6.0cm	3.5-6.2cm
	SRE-LR-200S	SRE-LR-200	SIE-LR-126S	SIE-LR-126	200-400cc	12.0cm	6.0cm	3.5-6.8cm
	SRE-LR-250S	SRE-LR-250	SIE-LR-127S	SIE-LR-127	250-500cc	12.0cm	7.0cm	3.8-7.4cm
	SRE-LR-340S	SRE-LR-340	SIE-LR-157S	SIE-LR-157	340-680cc	15.0cm	7.5cm	3.8-8.0cm
	BTE070	SPR070	n/a	n/a	70-90cc	6.0cm	4.0cm	3.5-4.6cm
	BTE250	SPR250			250-300cc	10.0cm	6.0cm	5.0-6.3cm
	BTE350	SPR350			350-420cc	12.0cm	6.0cm	5.0-7.0cm
	BTE500	SPR500			500-600cc	12.0cm	7.0cm	6.0-7.8cm
	BTE750	SPR750			750-900cc	15.0cm	7.5cm	7.0-9.0cm
	BTE1000	SPR1000			1000-1200cc	17.0cm	8.0cm	8.0-10.2cm
	SRE-CR-045S	SRE-CR-045	SIE-CR-045S	SIE-CR-045	45-90cc	8.7cm	5.5cm	2.7-4.5cm
	SRE-CR-100S	SRE-CR-100	SIE-CR-100S	SIE-CR-100	100-200cc	10.5cm	7.4cm	3.3-5.8cm
	SRE-CR-200S	SRE-CR-200	SIE-CR-200S	SIE-CR-200	200-400cc	13.7cm	10.5cm	4.0-7.0cm
	SRE-CR-350S	SRE-CR-350	SIE-CR-350S	SIE-CR-350	350-700cc	15.0cm	11.2cm	5.5-8.5cm

*No reinforced base on this device.


M172D 10/08



SPECIALTY SURGICAL PRODUCTS, INC.

Softspan™ Extremity Expander Wall Chart

Page 2 of 2



Expander Shape	Remote Port Cat #		Integrated Port Cat #		Suggested Fill Volume	A Length	B Width	C Projection
	Smooth Shell	Textured Shell	Smooth Shell	Textured Shell				
Round	SRE-RD-220S	SRE-RD-220	n/a	n/a	220-270cc	9.4cm	9.4cm	4.5-5.9cm
	SRE-RD-270S	SRE-RD-270			270-320cc	10.2cm	10.2cm	4.8-6.0cm
	SRE-RD-320S	SRE-RD-320			320-380cc	10.8cm	10.8cm	5.0-6.4cm
	SRE-RD-380S	SRE-RD-380			380-430cc	11.3cm	11.3cm	5.4-6.5cm
	SRE-RD-430S	SRE-RD-430			430-500cc	11.5cm	11.5cm	5.6-6.9cm
	SRE-RD-500S	SRE-RD-500			500-640cc	12.2cm	12.2cm	6.0-8.0cm
	SRE-RD-640S	SRE-RD-640			640-750cc	13.1cm	13.1cm	6.7-8.2cm
	SRE-RD-750S	SRE-RD-750			750-860cc	14.0cm	14.0cm	6.9-8.5cm
	SRE-RD-860S	SRE-RD-860			860-1020cc	14.5cm	14.5cm	7.7-9.0cm
	Curvilinear	SRE-CV-320S	SRE-CV-320	n/a	VU-CTE-320	320-480cc	19.0cm	6.5cm
SRE-CV-520S		SRE-CV-520		VU-CTE-520	520-780cc	23.8cm	7.5cm	4.5-7.0cm
SRE-CV-720S		SRE-CV-720		VU-CTE-720	720-1080cc	26.2cm	9.0cm	4.7-7.5cm
Angle	n/a	n/a	n/a	VU-ATE-200	200-300cc	19.0cm	4.0cm	3.5-4.5cm
				VU-ATE-450	450-650cc	23.5cm	5.0cm	4.5-6.0cm
Wedge	SRE-WG-100S*	SRE-WG-100*	n/a	n/a	100-120cc	8.3cm	5.5cm	4.5-5.5cm
	SRE-WG-250S*	SRE-WG-250*			250-300cc	10.5cm	7.5cm	6.0-7.3cm
Calf	SRE-CC-500S	n/a	n/a	n/a	500-600cc	19.0cm	8.0cm	5.5-7.0cm

*No reinforced base on this device.

ACCESSORIES

Catalog #	Description
DET01	Dermaspan™ Expander Templates
SET01	Softspan™ Expander Templates
PPL01	Precision Point™ Pendulum Locator
WNIS01	Winged Needle Infusion Set

Specialty Surgical Products, Inc. · Toll Free 888.878.0811

1131 North U.S. Hwy 93 · Victor, MT 59875 · email: info@ssp-inc.com

406.961.0102 · Fax: 406.961.0103