




SPECIALTY SURGICAL PRODUCTS, INC.

Dermaspan™ Low Pole Plus Expander Wall Chart



DERMASPAN™ LOW POLE PLUS (MID HEIGHT)

Textured Shell W/Lateral Extensions Catalog #	Volume	A* Width	B* Length	C Projection
LPP11	250-300cc	11.0cm	9.0cm	5.0-5.7cm
LPP12	310-370cc	12.0cm	9.8cm	5.2-6.0cm
LPP13	380-450cc	13.0cm	10.6cm	5.5-6.3cm
LPP14	460-550cc	14.0cm	11.2cm	6.0-7.3cm
LPP15	500-600cc	15.0cm	11.5cm	6.4-7.9cm
LPP16	580-690cc	16.0cm	12.0cm	6.6-8.2cm

DERMASPAN™ LOW POLE PLUS (FULL HEIGHT)

Textured Shell W/Lateral Extensions Catalog #	Volume	A* Width	B* Length	C Projection
LPP-FH11	250-300cc	11.0cm	10.2cm	5.0-5.7cm
LPP-FH12	310-370cc	12.0cm	11.0cm	5.2-6.0cm
LPP-FH13	360-430cc	13.0cm	11.9cm	5.4-6.3cm
LPP-FH14	480-575cc	14.0cm	12.9cm	6.1-7.3cm
LPP-FH15	600-720cc	15.0cm	13.8cm	6.4-7.9cm
LPP-FH16	750-900cc	16.0cm	14.7cm	6.9-8.5cm

*Nominal dimensions at minimum recommended fill volumes.